

Trouble-Free Training



ADAPTIVATION
The World as a Track


SwitchBoard

- Six switches in one confined area
- Each activation area is 4" x 5"
- 2 ounces of pressure to activate
- Can be used with any device that has 1/8th-inch jacks



Pal Pads

- Can be plugged into ANY communication aid or battery-operated assistive device
- Only need 1.2 ounces of force to activate
- 3 sizes
- Easily attach objects and textures
- Great for activities
 - Adapting books
 - Attaching to cases



Taction Pads

- Not a true switch
 - Has to be used with an Adaptation product
- Activated by moisture
- Clear, adhesive-backed switch
- Can be cut!



Flexible Switch

- Flat switch
- Adhesive-backed
- Pressure switch
- Works great in wheelchair headrests
- Flexible to the diameter of a coffee cup



Pneumatic Switch

- Activates with a puff of air
- Great way to make an object-based switch
- Can use squeeze toys, balls or bottles



7

Chipper

- Single-message device
- 20 seconds of record time
- External switch jack
- External speaker jack
- Relay jack for cause-and-effect



10

Lex

- Three levels
- One message per level
- 20 seconds of record time for each level
- 3 volume levels
- Pairs well with TMW Object Communication System



8



11

TalkingMyWay Object Communication System

- Encourages transitions from objects to symbols and/or text
- Great for schedules and choice making
- Child, Teen/Adult and Self-Care sets



9



12

Sequencer

- Plays messages in sequential order
- 2 levels (independent levels!)
- 45 seconds of record time per level (short play) or 90 seconds of record time per level (long play)
- External switch jack
- Relay jack for cause-and-effect



16



14



17



15



18



19



22

Sequencing Ideas

- Tasks/Directions/Steps
- Plays
- Singing Songs
- Jokes
- Introducing people
- Ordering at a restaurant
- Requesting assistance
- Checking out library books
- Steps to make photo copies
- Listening device/memory tool

20

Randomizing Ideas

- Choose reinforcer (good job, well done, great!)
- Greeting people (church, events, jobs)
- Assigning chores
- Conversation starters
- Games
- Cheering at a sporting event
- Choose your own adventure

23

Randomizer

- Plays messages in random order
- 2 levels (independent levels!)
- 45 seconds of record time per level (short play) or 90 seconds of record time per level (long play)
- External switch jack
- Relay jack for cause-and-effect



21

VoicePal 8

- Up to 8 messages can be recorded and played using adaptive switches
- Delayed Activation
- 90 seconds record time (11 per message)



24



25

VoicePal 8K

- Up to 8 messages can be recorded or played using adaptive switches or the on-board keypad
- Delayed Activation
- Repeat Message Mode
- 90 seconds of record time (11 seconds each)



28



26



29



27



30

VoicePal Levels

- 4 Levels
- 17.5 minutes of record time
- 10 external switch jacks
- On-board keypad
- Configure for 2, 3, 5 or 10 messages
- Scanning capabilities
- Sequencing mode
- Randomizing mode



Battery Interrupter



LinkSwitch

- Used to activate battery-operated devices using switches
- Dual-channel modes
- Modes of control: Direct, Latch, Timing (seconds or minutes) and Alarm



FreeHand

- Wireless control of AC appliances (using X-10 module)
- Modes of control: Latch, Direct and Timing (seconds or minutes)
- Range: 75 feet



37


Resources

- Recipes for Success Remixed
 - Sample recipes available at www.adaptivation.com
- Application Photo Album of Ideas
- Workshop Kit
 - (includes shipping and insurance both ways in the Continental US)
- Trial Program
 - (includes shipping and insurance one way in the Continental US)

40



38



39